

Levo ME Pipette Filler

Levo ME Pipette Filler is a new design that has better ergonomics in pipetting and can fit almost all commonly used glass or plastic pipettes.

01-2210

Features

- Lightweight electronic pipette filler, easy for one-hand operation
- Fits almost all commonly used glass or plastic pipettes, from the volume of 0.1-100mL
- Easy to adjust the variable speed
- Hydrophobic autoclavable filter prevents over-aspiration
- Stable and reliable performance, elegant appearance



Specifications

	2210
Aspirate Speeds	continuous
Dispense Speeds	continuous
Battery	Lithium-Ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2~3 Hours
Pipette Types	Glass or Plastic Pipette(0.1-100mL), Pasteur Pipettes
Filter	0.45 µm Hydrophobic